REGO-LOK™ for Securing CGA Fittings on Liquid Cylinders

Application

The REGO-LOK[™] is designed for installation on the RegO T9450 and T9460 Series liquid cylinder valves to deter and prevent the removal of the CGA fitting from the valve. The REGO-LOK[™] retains standard CGA outlet connection so unauthorized persons do not remove the fitting. By use of a special one-way bolt, the REGO-LOK[™] is secured to the valve. The REGO-LOK[™] installs in a few minutes with the use of screwdrivers, without valve disassembly, brazing, welding, or drilling. The REGO-LOK[™] deters and prevents fitting removal by gas customers, however allows the replacement of fittings by authorized gas supplier plant personnel.

Use The REGO-LOK[™] for compliance with CGA SB-26 for medical and industrial liquid cylinders.

Features

- Stainless Steel REGO-LOK[™] with one-way bolt
- · Retrofit all common liquid cylinder valves
- · Can be supplied on new RegO liquid cylinder valves
- REGO-LOK™ indicates "WARNING: DO NOT REMOVE"
- Worn CGA fittings can be simply replaced by authorized personnel. Requires new 9464RL-6 Bolt
- Can fit over existing fittings for CGA 540, CGA 440, CGA 295, CGA 320, and CGA 326. Check fitting hex size

NOTE: RegO supplied fitting P/N CGA580RL is required for REGO-LOKTM use with CGA 580 connection


· Prevents loosening of CGA fittings on valves





REGO

RegO-Lok™

Satisfies CGA SB-26 and FDA requirements for medical and industrial liquid cylinders.







Ordering Information

Part Number	Item Description	Typical Service Connection
9464RL-0	REGO-LOK™ for ¾" hex fittings	N/A
9464RL-1	REGO-LOK™ for ⁷ ⁄₀" fittings	CGA 320, CGA 326 & CGA 295
9464RL-2	REGO-LOK™ for 1" fittings	CGA 440, CGA 540
9464RL-3	REGO-LOK™ for 11⁄8" hex CGA 580RL fitting by RegO	CGA 580
CGA580RL	3/6" MNPTxCGA for use with 9464RL-3	CGA 580

DO NOT REMOVE

